

Saving to

My Library

All About HRV Part 5: Data Editing – MindWare ...
Snapshot

Training Videos

Training Events

MindWare Training

Resources

Manuals

Downloads

Create a Support Case

Navigate

MindWare Website

Online Store

Resources

Search Support



All About HRV Part 5: Data Editing

Posted 26 September 2017 | By [Eric Morgan](#) | Under [Blog](#)

Select analysis version to view the applicable content: 3.2 ▼

We have now gone through the [various statistics](#) in the [HRV](#) application, how they are calculated, and [what they infer](#). To ensure that these statistics are accurate, the R peaks on the ECG signal need to be checked for correctness.

Effects of Artifact on RSA

It is critical that any and all artifact be corrected to ensure valid HRV-related statistics. This is true even more so than when [editing ECG for Impedance Cardiography](#), since a single erroneous peak can significantly change the results. Let's take for example this segment which is free from artifact:

Contents ▼

1 Effects of Artifact on RSA


▼ 2 Editing Technique

2.1 30 Second Rule

2.2 10% Rule

3 Edit Files

4 Wrapping Up



Make note of the RSA value. Now let's delete a single R peak



You will notice the RSA value increased by over 3 units! This is a statistically significant change in RSA, resulting from only a single missed beat. In this particular case, we can see the drastic spike in the IBI series which is not physiologically plausible, and need to insert a beat here where there should be one.

Editing Techniques


ECG editing has been gone over thoroughly in [previous blog posts](#). I won't reiterate everything here, but these are the unique requirements when editing ECG for HRV:

Saving to My Library ☐


All About HRV Part 5: Data Editing – MindWare ...
Snapshot

30 Second Rule

A minimum of 30 seconds of contiguous R peak data is needed to effectively calculate RSA. To start, be sure that you are looking at segments of time that are greater than 30 seconds. 60 second segments are recommended as a minimum in practice, as they offer the most flexible editing options. Remember that you can remove noisy data from the beginning or end of a segment so long as there is at least 30 seconds of data remaining in the middle



but it is **never** acceptable to delete a portion of data in the middle



This gap in the middle of the segment will be interpreted as a long IBI and will skew the calculations. This can be seen reflected in the IBI series



10% Rule


The **10% rule** is in reference to how many estimated beats are acceptable before the statistics are not longer valid. Beats are typically estimated during periods of signal noise in which an R peak cannot be identified, or when dealing with cardiac arrhythmia. The rule states that no more than 10% of total beats within a segment, epoch, or subject should be

estimated. If this limit is exceeded, it is best to discard that segment/epoch/subject.

Be sure to use the Point Distribution editing statistics in the hand corner of the analysis screen to keep track of how many edits have been made in the current segment.


Saving to My Library ☐

All About HRV Part 5: Data Editing – MindWare ...
Snapshot



Edit Files

The edit file stored by the HRV application has a **.edh2** file extension. They are typically saved in an **edit data** folder in the same directory as the data file.



It is also possible to [load an edit file](#) from the [Cardiac Impedance Analysis application](#), but if you do so make sure that the [edits are valid for HRV](#).

Wrapping Up

Now we have discussed the origins of HRV, the signals and processes needed to derive meaningful information in this domain, and techniques for ensuring our results are accurate and correct. With this information, you should be able to confidently and effectively use heart rate variability in your study.







Saving to






My Library


☐

All About HRV Part 5: Data Editing – MindWare ...
Snapshot


Tagged: # 10% Rule # 30 Second Rule # AllAboutHRV # ECG Editing
Editing


-  Support Portal Home
-  Knowledge Base
-  MindWare Blog
-  Manuals
-  Downloads
-  Create a Support Case


-  Training Guides
-  Premium Training
-  Training Videos
-  Training Events
-  QuickSupport Agent



Subscribe to our mailing list to keep up to date on news, webinars, and software updates!

 MindWare Homepage

 Web Store

 Privacy Policy

Submit Feedback
© Copyright MindWare Technologies, LTD.