

Kubios HRV

Heart Rate Variability Analysis Software

Mika Tarvainen, PhD
CEO, Kubios Oy




MATLAB EXPO 2018

15.5.2018 Helsinki

What is HRV

Regulation of heart rate




Autonomic nervous system (ANS)

- To preserve blood pressure homeostasis, ANS regulates
 - Heart rate and cardiac function
 - Vasomotor activity
 - Arterial baroreflex
- ANS divided into sympathetic and parasympathetic branches

HRV time series

- Continuous regulation of heart rate (HR) → Heart rate variability (HRV)
- Sympathetic activity (**GAS**)
→ Increases HR and decreases HRV
- Parasympathetic activity (**BRAKE**)
→ Decreases HR and increases HRV
- Respiratory sinus arrhythmia (RSA)
- High HRV indicates good recovery and high ability to tolerate stress







HRV Analysis Methods

- Time-domain (Mean RR, SDNN, RMSSD etc.)
 - Intensity of HRV
- Frequency-domain (LF and HF components, LF/HF ratio etc.)
 - Frequency components of HRV (assessment of sympathovagal balance)
- Nonlinear (Entropy measures etc.)
 - Complexity of HRV
- Time-varying
 - HRV dynamics

$$SDNN = \sqrt{\frac{1}{N-1} \sum_{j=1}^N (RR_j - \overline{RR})^2}$$

$$RMSSD = \sqrt{\frac{1}{N-1} \sum_{j=1}^{N-1} (RR_{j+1} - RR_j)^2}$$



Applications of HRV

Medical Research



- Risk evaluation after MI
- Evaluation of heart failure
- Marker of diabetic neuropathy
- Evaluation of CAD
- ANS testing
- Sleep apnea
- Affective disorders
- Anesthesia monitoring etc.

Stress and Wellbeing



- Occupational stress
- HRV Biofeedback
- Resonant frequency
- Physiotherapy, yoga etc.

Fitness and Exercise



- Fitness assessment
→ maxHR, anaerobic threshold, energy consumption
- Recovery monitoring
→ Avoiding overtraining
- Coaching
→ HRV based individualized training prescription

Kubios Oy

About the company

- Founded 2016
- Software and algorithms development
- Market leader in HRV analysis software for scientific research and professional use
- Customers: Researchers, wellbeing therapists, sports/exercise coaches and athletes, personal monitoring etc.



Mika Tarvainen, PhD
CEO

15 years of experience in medical signal analysis and physiological modeling



Jukka Lipponen, PhD
COO

Expert on ECG signal analysis and software development



Juha-Pekka Niskanen
CPO

Qualified Medical Physicist, Expert on Matlab software development




Perttu Ranta-aho
CTO

Expert on signal analysis, IT solutions and software development

Facts and figures

- Golden Standard HRV software
- Used in 1000 scientific studies
- Device independent software
- Over 40 analysis parameters
- Used in 120 countries (6/7 continents)



Products

Kubios HRV Standard (ver. 3.1)

- For non-commercial personal use
- Supports most common HR monitors (RR data)
- Standard HRV analysis
- Freeware




Kubios HRV Premium (ver. 3.1)

- For research and professional use
- Supports several HR monitors, ECG and PPG devices
- Full-featured HRV analysis software
 1. Extended data support (ECG and PPG data)
 2. Built-in beat detection
 3. Automatic correction of missed, extra and misaligned (ectopic) beats
 4. ECG derived respiration (EDR)
 5. Extended analysis features (spectrogram with “fire” colormap)
 6. Extended reports and exporting options (CSV batch file export)



HRV Premium


Kubios HRV Premium is the market leader in heart rate variability (HRV) analysis software for scientific research and professional use.



How we use MATLAB

How it all started?

- MATLAB 4.2 and 5 (1994-1999)
 - MSc studies in Medical Physics
 - Used in many courses for practical works and exercises
 - MATLAB 6, 7 and 8 (2000-2014)
 - PhD research in Medical Signal Analysis
 - Development of HRV analysis software started
- Ver. 1.1 released in Sep 2002



Ver. 1.1 (9/2002)

- ✓ GUI
- ✓ Standard HRV Analysis
- ✓ RR data support
- ✓ Windows

Ver. 2.0 (10/2008)


- ✓ Nonlinear HRV parameters
- ✓ Support for Polar and Suunto IBI data
- ✓ Windows & Linux

Ver. 2.1 (7/2012)

- ✓ QRS detector
- ✓ ECG data support
- ✓ Windows & Linux

Ver. 2.2 (5/2014)

- ✓ ECG derived respiration
- ✓ Updates to data support
- ✓ Windows, Linux & Mac



MATLAB Compiler

Limited functionality
(no MCR)

MATLAB Compiler Runtime

Still relying on MATLAB

1. Optimal for algorithm development and testing

Example: Signal baseline fitting

$$\hat{\theta}_{\lambda} = \arg \min_{\theta} \{ \|H\theta - z\|^2 + \lambda^2 \|D_d(H\theta)\|^2 \}$$


$$\hat{\theta}_{\lambda} = (H^T H + \lambda^2 H^T D_d^T D_d H)^{-1} H^T z$$

$$\hat{z}_{\text{trend}} = H \hat{\theta}_{\lambda}$$

$$\hat{z}_{\text{stat}} = z - H \hat{\theta}_{\lambda} = (I - (I + \lambda^2 D_2^T D_2)^{-1}) z$$

$$H = I \in \mathbb{R}^{(N-1) \times (N-1)}$$

$$D_2 = \begin{pmatrix} 1 & -2 & 1 & 0 & \cdots & 0 \\ 0 & 1 & -2 & 1 & \ddots & \vdots \\ \vdots & \ddots & \ddots & \ddots & \ddots & 0 \\ 0 & \cdots & 0 & 1 & -2 & 1 \end{pmatrix} \in \mathbb{R}^{(N-3) \times (N-1)}$$



2. Graphics

MATLAB figure

- ✓ uimenu and toolbar
- ✓ 2D & 3D visualisations
- ✓ Highly customizable controls for GUI design
- ✓ Mouse/keyboard interactions

MATLAB UI figure (App designer)

- ✓ Interesting, but still some limitations (mouse/keyboard interaction, UImenu support)



3. Fast to build applications

MATLAB Compiler

- ✓ Matlab code → Windows, macOS and Linux applications
- ✓ deploytool – easy to use, but we need to customize our installers
- ✓ Code signing
- ✓ Agile software development



Kubios HRV Premium – DEMO

✓ Export data from your measurement device



Kubios HRV Premium – DEMO


✓ Export data from your measurement device



Suunto HR monitors (www.suunto.com)

- ✓ T6 series and Memorybelt
- ✓ Ambit 1, 2 and 3 series
- ✓ Spartan series

→ Export FIT file from Movescout



Kubios HRV Premium – DEMO

- ✓ Export data from your measurement device



Garmin HR monitors (www.garmin.com)

- ✓ Forerunner 235, 620, 630, 920XT, 735 XT, 935
- ✓ Fenix 3, 3HR, and 5; Edge 520, 820, 1000

NOTE: Turn on HRV recording by changing “Log HRV” setting (Fenix 5 and Forerunner 935) or follow the steps give in <https://sporttracks.mobi/blog/tracking-hrv-garmin-watches>

→ Export FIT file from Garmin Connect



Kubios HRV Premium – DEMO

✓ Export data from your measurement device



Polar HR monitors (www.polar.com)

✓ V800, RS800, RS800CX, CS600, S810

→ Perform RR recording test and export RR data from Polar Flow web service




Kubios HRV Premium – DEMO

- ✓ Export data from your measurement device

Supported ECG and PPG devices

- ✓ Actiheart and Actiwave Cardio (www.camntech.com)
- ✓ Alivecor Kardia (www.alivecor.com)
- ✓ Biopac ECG and PPG devices (www.biopac.com)
- ✓ Bittium Faros ECG (www.bittium.com)
- ✓ Empatica E4 (www.empatica.com)
- ✓ FirstBeat Bodyguard (www.firstbeat.com)
- ✓ Mindfield MindMaster (www.mindfield.de)
- ✓ Shimmer ECG and PPG devices (www.shimmersensing.com)
- ✓ Zephyr Bioharness (www.zephyranywhere.com)


etc.



Kubios HRV Premium – DEMO

✓ Analysing data in Kubios HRV

1) Opening a recording



Kubios HRV Premium – DEMO

✓ Analysing data in Kubios HRV

- 1) Opening a recording
- 2) Check beat detection and correct artefacts if necessary



Kubios HRV Premium – DEMO

✓ Analysing data in Kubios HRV

- 1) Opening a recording
- 2) Check beat detection and correct artefacts if necessary
- 3) Place as many analysis samples as you want (select stationary time periods)
- 4) All HRV analysis results are computed and visualised immediately



Kubios HRV Premium – DEMO

✓ Analysing data in Kubios HRV


- 1) Opening a recording
- 2) Check beat detection and correct artefacts if necessary
- 3) Place as many analysis samples as you want (select stationary time periods)
- 4) All HRV analysis results are computed and visualised immediately
- 5) Apply time-varying analysis



Kubios HRV Premium – DEMO

✓ Analysing data in Kubios HRV


- 1) Opening a recording
- 2) Check beat detection and correct artefacts if necessary
- 3) Place as many analysis samples as you want (select stationary time periods)
- 4) All HRV analysis results are computed and visualised immediately
- 5) Apply time-varying analysis
- 6) Save results
→ PDF, TXT, MAT + “SPSS friendly” batch file



Kubios HRV Premium – DEMO

✓ Kubios HRV Premium reports


1) Report pages (1 page/sample)



Kubios HRV Premium – DEMO


✓ Kubios HRV Premium reports

- 1) Report pages (1 page/sample)
- 2) Report page for time-varying analysis
- 3) CSV-file
- 4) MAT-file
- 5) “SPSS friendly” batch file



Thank you

Mika Tarvainen
CEO, Kubios Oy
mika.tarvainen@kubios.com
www.kubios.com



Follow us:   